



Bulb syringes

with mat surface



Size	Volume in cm ³ (ml)	Opening in mm	Item no.*	D in mm	L in mm	D1 in mm
0	8	6	4100.2.xx000	30	39	6
2	25	6	4102.2.xx000	42	55	6
3	40	6	4103.2.xx000	47	63	6
5	89	6	4105.2.xx000	62	79	6
7	143	6	4107.2.xx000	71	95	6
9	224	6	4109.2.xx000	83	111	6
11	347	6	4111.2.xx000	93	133	6
13	483	6	4113.2.xx000	103	150	6
9	224	14	4119.2.xx000	83	111	14
11	347	14	4121.2.xx000	93	133	14
13	483	16	4133.2.xx000	103	150	16

* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

Standard colors:

00 red-brown, 10 pastel blue and 20 black,
other colors, special openings and special weights upon request.

Bulb material:

SAROLIT™
physiologically unobjectionable, phthalate-free, cadmium-free,
latex-free, Shore hardness 52 A, other Shore hardness upon request.