


Bulbs with adapters

different combination options




5103.2.10165



7203.2.10265



7323.2.10788



Name	Volume in cm ³ (ml)	Item no.*	Adapter	L in mm	D in mm
Inline pump size 3 with inlet and outlet valve (red/green)	65	5103.2.xx165	Z632.03 + Z632.40	93	49
Pump bulb with tubing connector for caulyt bulbs	45	7203.2.xx265	Z634.20	82	43
Pump bulb standard with tubing connector	48	7323.2.xx788	Z241.30	88	41

* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

All adapter can be combined with different pump and suction bulbs.
Our sales team would be glad to advise you.

Standard colors: bulbs in 00 red-brown, 10 pastel blue, 20 black,
other colors upon request.

Bulb material: SAROLIT™
physiologically unobjectionable, phthalate-free, cadmium-free, latex-free,
Shore hardness 52 A, other Shore hardness upon request.