


Bulbs with pressure control valves

with mat surface




5103.2.10142



7203.2.10151



7323.2.10153



Name	Volume in cm ³ (ml)	Item no.*	Pressure control valve	Material	L in mm	D in mm
Inline pump size 3 with pressure control valve Z652.20	65	5103.2.xx142	Z652.20	Plastic with metal insert	93	49
Pump bulb with pressure control valve Z674.20	45	7203.2.xx151	Z674.20	Metal (chrome-plated)	82	43
Pump bulb standard with pressure control valve Z640.20	48	7323.2.xx142	Z640.20	Plastic	88	41

* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

All pressure control valves can be combined with different pump and suction bulbs.
Our sales team would be glad to advise you.

Standard colors: bulbs in 00 red-brown, 10 pastel blue, 20 black,
other colors upon request.

Bulb material: SAROLIT™
physiologically unobjectionable, phthalate-free, cadmium-free, latex-free,
Shore hardness 52 A, other Shore hardness upon request.