



Pump bulb flattened



Name	Volume in cm ³ (ml)	Item no.*	L in mm	D in mm
Pump bulb flattened with one inserted valve	28.5	1073.2.xx001	91	43

* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

Standard colors: 00 red-brown, 10 pastel blue and 20 black, other colors upon request.

Bulb material: SAROLIT™
physiologically unobjectionable, phthalate-free, cadmium-free,
latex-free, Shore hardness 52 A, other Shore hardness upon request.