


Douche syringes standard size 13



1-part set

2-part set

Name	Head type	Volume in cm ³ (ml)	Item no.*	L in mm	D in mm	L1 in mm
Vaginal douche size 13, mat, 1-part set	Standard**	483	4133.2.xx080	150	103	270
Vaginal douche size 13, mat, 1-part set	Rotation***	483	4133.2.xx081	150	103	270
Vaginal douche size 13, mat, 2-part set	Standard**	483	4133.2.xx082	150	103	273
Vaginal douche size 13, mat, 2-part set	Rotation***	483	4133.2.xx083	150	103	273
Vaginal douche size 13, matt, pipe for douche syringes	Rotation***	483	4133.2.xx277	150	103	264

* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

** Standard head: Liquid emerges through 6 openings at the top. *** Rotation head: Liquid emerges rotating at the top.

Standard colors:

00 red-brown, 10 pastel blue, 20 black, special colors upon request.

Bulb material:

SAROLIT™
physiologically unobjectionable, phthalate-free, cadmium-free
Shore hardness 52 A, other Shore hardness upon request.

The vaginal douches are standardly delivered with collar (splash guard) made of SAROLIT™. Also available without collar upon request. The vaginal douche sets are made of polyethylene (PE).
1-part set means – only vaginal douche collar (Ø approx. 19 mm),
2-part set means – a rectal pipe is additionally located below the vaginal douche collar, which can be removed.

