


# Ear syringes

with mat surface



Size	Volume in cm <sup>3</sup> (ml)	Item no.*	L in mm	D in mm
0	8	8100.2.xx000	58	29,5
3	27	8103.2.xx000	84	44
5	41	8105.2.xx000	101	50
7	75	8107.2.xx000	117	58
9	135	8109.2.xx000	138	71
11	206	8111.2.xx000	148	83
13	317	8113.2.xx000	173	93

\* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

**Standard color:**

00 red-brown, 10 pastel blue, 20 black,  
other colors and special weights upon request.

**Bulb material:**

SAROLIT™  
physiologically unobjectionable, phthalate-free, cadmium-free,  
latex-free, Shore hardness 52 A, other Shore hardness  
upon request.

