



# Ear syringes

US sizes with mat surface



Size	Volume in cm <sup>3</sup> (ml)	Item no.*	L in mm	D in mm
1-oz	33	8201.2.xx000	90	44

\* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

**Standard color:**

00 red-brown, 10 pastel blue, 20 black,  
other colors and special weights upon request.

**Bulb material:**

SAROLIT™  
physiologically unobjectionable, phthalate-free, cadmium-free, latex-free,  
Shore hardness 52 A, other Shore hardness upon request.