




# Nasal aspirator bulbs



P200.0.10000



P250.0.10000



Name	Volume in cm <sup>3</sup> (ml)	Material	Item no.*	L in mm	D in mm
Nasal aspirator bulb (neutral)	34	SAROLIT™	P200.0.xx000	59	41
Nasal aspirator bulb (neutral)	34	TPE-S	P400.0.xx000	59	41
Nasal aspirator bulb with rings	34	SAROLIT™	P250.0.xx000	59	42
Nasal aspirator bulb with rings	34	TPE-S	P450.0.xx000	59	42

\* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

**Standard colors:** 10 pastel blue and 30 white,  
other colors upon request.

**Bulb material:** SAROLIT™  
physiologically unobjectionable, phthalate-free, cadmium-free,  
latex-free, Shore hardness 52 A, other Shore hardness upon request.

TPE-S  
physiologically unobjectionable, phthalate-free, cadmium-free,  
latex-free, Shore hardness 60 A, other Shore hardness upon request.