


Nasal aspirators

with nasal olive



Name	Volume in cm ³ (ml)	Item no.*	L in mm	D in mm
Pump bulb with nasal olive	40	7203.2.xx053	100.5	43

* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

Standard colors: bulb in 00 red-brown, 10 pastel blue, 20 black, other colors upon request.

Bulb material: SAROLIT™
physiologically unobjectionable, phthalate-free, cadmium-free, latex-free,
Shore hardness 52 A, other Shore hardness upon request.