



Pump bulbs



Name	Volume in cm ³ (ml)	Item no.*	L in mm	D in mm	D1 in mm
Standard with one inserted valve	47	7323.2.xx001	89	41	6
Standard with one inserted valve	47	7313.2.xx001	89	41	7

* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

Standard colors: 00 red-brown, 10 pastel blue and 20 black, other colors upon request.

Bulb material: SAROLIT™
physiologically unobjectionable, phthalate-free, cadmium-free, latex-free, Shore hardness 52 A, other Shore hardness upon request.