


# Travel douches

with flexible pipe



Name	Volume in cm <sup>3</sup> (ml)	Pipe material	Item no.*	L in mm	D in mm	L1 in mm
Vaginal douche bulb 10-oz with flexible white pipe	330	SAROLIT™	4010.3.xx225	136	78	264.5
Vaginal douche bulb 10-oz with flexible white pipe and collar	330	SAROLIT™	4010.3.xx223	136	78	264.5

\* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

**Standard colors:** bulb in 00 red-brown, 10 pastel blue, 30 white, vaginal douche pipe 20 black, 30 white, other colors upon request.

**Bulb material:** SAROLIT™  
physiologically unobjectionable, phthalate-free, cadmium-free, latex-free, Shore hardness 52 A, other Shore hardness upon request.

For the above mentioned vaginal douches, the vaginal douche pipe can be inserted or screwed into the bulb for space-saving reasons. The dimensions in that case are: For item 4010... instead of 270 x Ø 81 mm -> 160 x Ø 81 mm.