


Douche syringes standard size 9



1-part set

2-part set

Name	Head type	Volume in cm ³ (ml)	Item no.*	L in mm	D in mm	L1 in mm
Vaginal douche size 9, mat, 1-part set	Standard**	224	4119.2.xx080	111	83	231
Vaginal douche size 9, mat, 1-part set	Rotation***	224	4119.2.xx081	111	83	231
Vaginal douche size 9, mat, 2-part set	Standard**	224	4119.2.xx082	111	83	234
Vaginal douche size 9, mat, 2-part set	Rotation***	224	4119.2.xx083	111	83	234

* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

** Standard head: Liquid emerges through 6 openings at the top. *** Rotation head: Liquid emerges rotating at the top.

Standard colors: 00 red-brown, 10 pastel blue, 20 black, special colors upon request.

Bulb material: SAROLIT™
physiologically unobjectionable, phthalate-free, cadmium-free
Shore hardness 52 A, other Shore hardness upon request.

The vaginal douches are standardly delivered with collar (splash guard) made of SAROLIT™. Also available without collar upon request.

The vaginal douche sets are made of polyethylene (PE).

1-part set means – only vaginal douche collar (Ø approx. 19 mm),

2-part set means – a rectal pipe is additionally located below the vaginal douche collar, which can be removed.

