


# Pipette bulb



Name	Volume in cm <sup>3</sup> (ml)	Item no.*	L in mm	D in mm	D1 in mm
Pipette bulb	5	3100.2.xx000	39	24	5

\* The "xx" in the item numbers is replaced with the color codes (see color chart on fold-out).

**Standard colors:** 00 red-brown, 10 pastel blue and 20 black,  
other colors, special openings and special weights upon request.

**Bulb material:** SAROLIT™  
physiologically unobjectionable, phthalate-free, cadmium-free, latex-free,  
Shore hardness 52 A, other Shore hardness upon request.

**Bulb surface:** mat